

## Compact Shaker KS 15 B



### Description

The reliable Horizontal Shaker KS 15 B for general shaking tasks for daily use. The loading capacity is 15 kg on a shaking platform of 400 x 300 mm.

### Basic Equipment

Basic device without shaking plate, **without rack system**

### Loading Capacity

Erlenmeyer flasks 100 ml	20 pieces
Erlenmeyer flasks 250 ml	12 pieces
Erlenmeyer flasks 500 ml	6 pieces
Erlenmeyer flasks 1000 ml	4 pieces
Separating funnels 250 ml	1 pieces

You will find further loading possibilities [here](#)

## Details

<b>Motion:</b>	reciprocating (horizontal)
<b>Shaking platform:</b>	400 x 300 mm
<b>Max. load:</b>	15 kg
<b>Shaking speed:</b>	30 - 420 rpm, steplessly adjustable
<b>Stroke:</b>	17 mm
<b>Timer:</b>	0 - 120 min / continuous
<b>Order No.:</b>	6173 000
<b>Accessory:</b>	Rack System Combifix KS Universal tray KS Incubator Hood TH 15

## Technical specifications

<b>Electrical supply:</b>	230 V or 115 V, 50 / 60 H please indicate in case of order
<b>Enclosure protection:</b>	IP 21
<b>Heat emission:</b>	approx. 5 - 10 W
<b>Ambient temperature:</b>	5°C to 50°C
<b>Relative humidity:</b>	~ 85 %
<b>Dimensions (w x d x h):</b>	510 x 490 x 150 mm
<b>Weight:</b>	18 kg