

# Three-Dimensional Shaker TL 2



## Description

The Laboratory Shaker TL 2 works with slow, three-dimensional motion for gentle, uniform wetting. Well suited for petri dishes, culture bottles etc.

### Basic Equipment

Basic device incl. basic tray and rubber mat

## Details

- Motion:** three-dimensional orbital motion  
**Shaking platform:** 220 x 220 mm  
**Max. load:** 2 kg

**Shaking speed:** 1 - 22 rpm, steplessly adjustable  
**Angle of tilt:** 8° fixed  
**runtime:** 0 - 9 h 59 min / continuous  
**Order No.:** 6139 000  
**Accessory:** Additional Tray TL 2

## Technical specifications

**Electrical supply:** 100 - 240 V, 50 / 60 Hz  
**Enclosure protection:** IP 21  
**Heat emission:** max. 3 W  
**Ambient temperature::** 5°C to 45°C  
**Relative humidity:** ~ 85 %  
**Dimensions (w x d x h):** 250 x 415 x 118 mm  
**Weight:** 9,1 kg