

## Tilt Shaker WS 10



### Description

The Tilt Shaker WS 10 works with slow, linear tilting motion and a fix angle of tilt of 5° with a shaking speed of 1 - 23 rpm. The shaking platform of 340 x 340 mm can be loaded with a weight of 10 kg.

### Basic Equipment

Basic device incl. basic tray and rubber mat

### Details

**Motion:** linear tilting motion

**Shaking platform:** 340 x 340 mm

**Max. load:** 10 kg

**Shaking speed:** 1 - 23 rpm, steplessly adjustable  
**Angle of tilt:** 5° fixed  
**Timer:** 0 - 120 min / continuous  
**Order No.:** 6114 000  
**Accessory:** Multi-storey rack system WS / TL

## Technical specifications

**Electrical supply:** 24 V DC, 50 / 60 Hz 115 V version on request  
**Enclosure protection:** IP 21  
**Heat emission:** approx. 8 W  
**Ambient temperature:** 5°C to 50°C  
**Relative humidity:** ~ 85 %  
**Dimensions (w x d x h):** 355 x 450 x 205 mm  
**Weight:** 18,5 kg