



ChemFAST Classic ChemFAST Sharp

Constant air volume and variable air volume. Total exhaust fume cupboards



TOTAL EXHAUST TECHNOLOGY	EN 14175	CONSTANT/ VARIABLE AIR VOLUME	WIDTH (m) 1.2 1.5 1.8 2.1 2.4
--------------------------	----------	-------------------------------------	--

ChemFAST Classic are total exhaust chemical fume cupboards working at **Constant Air Volume technology (CAV)**. The inflow velocity can be adjusted from 0,3 to 0,7 m/sec, both during the installation and the ordinary maintenance; this feature makes it possible to choose the safest velocities for the removal of the contaminants.

ChemFAST Sharp are total exhaust chemical fume cupboards working at **Variable Air Volume technology (VAV)**, which regulates the volume of exhausted air depending upon the opening of the sliding front door.

Both systems assure flexibility, safety and user friendliness and represent the ideal solution for chemical labs thanks to a combination of high standards of protection for the operator with a reduced air consumption.